High throughput orientation

mPrep/s™ Workstation


orient specimens quickly and precisely in mPrep/s capsules — without embedding molds or further handling






mPrep/s Workstation

The mPrep/s Workstation™ is designed for efficient workflow for mounting biological and material specimens. It has areas for holding mPrep/s capsules and their removable screens, a large flat working surface for specimen prepar-







The removable screens are easy to load onto the insertion tool in the center of the Workstation. It holds them horizontally, simplifying placement of specimens in the desired orientation for processing in mPrep/s capsules. Specimens remain oriented throughout processing and embedding in capsules. End-on mounting is achieved by opening a slit in the top screen using the Workstation actuator lever, gently clamping specimens in place vertically. This is useful for mounting fibers, nerves, muscle, plant tissues and other specimens for subsequent cross-sectioning or for vertical mounting for SEM.

Features of the mPrep/s Workstation™

- Layout enhances workflow efficiency and prevents errors.
- Dissection wells near the center for preparing specimens while keeping them moist. Sloped sides provide easy clean out and prevent specimens from being trapped in corners.
- High density polyethylene (HDPE) work surface is easy to clean. Minimizes dulling of fine dissection tools. Provides low glare/high contrast viewing of specimens. Top is replaceable if it becomes marred with cuts after long use.
- Entire Workstation made from materials for long service and durability — even in contact with saline and buffer solutions.
- Rounded front edge provides comfortable hand rest.

