


Automated Immunogold Labeling of TEM Grids

mPrepTM
System

Immunogold labeling (IGL) was introduced in the 1970s to enable electron microscopists to localize the constituents of a cell by attaching (electron dense) antibody-bound colloidal gold particles to antigens on specific cellular components. Here we demonstrate use of the Microscopy Innovations mPrepTM System ASP-1000 Automated Specimen Processor for this application.


Scientists at the University of Wisconsin Medical School Electron Microscope Facility used the mPrepTM System ASP-1000 Automated Specimen Processor to perform unattended, fully-automated IGL grid preparation using an AurionTM blocking buffer and protocol with excellent results.



ASP-1000 Automated Specimen Processor

Thin sections of resin-embedded, nematode samples (*C. elegans*) were placed on TEM grids, which were then inserted into mPrep/gTM capsules. The capsules were mounted on the ASP-1000 and processed unattended to accomplish IGL in 7.5 hours. The following Aurion protocol was utilized:


Reagent	Incubation Time (minutes)	Rinsing Time (minutes)	Number of Rinsing Steps
Glycine	15	n/a	n/a
Aurion Blocking Solution	15	5	3
Primary Antibody	120	5	6
Secondary Antibody (gold conjugate)	120	5	9
Glutaraldehyde	5	5	9



MICROSCOPYTM
INNOVATIONS LLC

213 Air Park Rd Suite 101
Marshfield WI 54449-8626
888.302.3925 or 715.384.3292
info@microscopyinnovations.com


Technical Note

**For more information, contact:**

Thomas E. Strader, MS-Biotech
Automation Engineer
608.770.7649
tom.strader@microscopyinnovations.com

Reagents on the ASP-1000 deck were stored in 96-well sealed microplates.

Rinsing reagents were pre-dispensed into 96-well polypropylene microplates and sealed with aluminum film. Glycine, glutaraldehyde and antibodies (primary and secondary gold conjugate) were added to the microplate immediately prior to starting the pre-programmed protocol. Samples were automatically agitated by dispensing and immediately re-aspirating 40 μ L every 10 minutes during antibody incubations and every 15 seconds during rinsing steps. Upon protocol completion, transmission electron microscopy (TEM) was performed at the UW Medical School EM Facility.



A component of the endosomal sorting complex required for transport (ESCRT) machinery localizes to the limiting membrane of a multivesicular endosome in a *C. elegans* 1-cell stage embryo, as revealed by immunogold labeling (6 nm Au particles). Scale bars = 200 nm.

Acknowledgments: Special thanks to Ben August (University of Wisconsin - Madison School of Medicine and Public Health EM Facility) for assistance in adapting [the Aurion post-embedding immunogold labeling protocol](#) to an automated version and to Elisa Frankel of the Audhya lab in the UW - Madison Department of Biomolecular Chemistry for providing specimens to process and image.

Rev 1 160217