



## Quick Start — BSL3/4 Virus Preparation

### Biosafety level 3/4 containment suite



### Electron microscopy laboratory



**MICROSCOPY**  
INNOVATIONS<sup>LLC</sup>

213 Air Park Rd Suite 101  
Marshfield WI 54449-8626  
888.302.3925 or 715.384.3292  
info@microscopyinnovations.com



## Quick Start — BSL3/4 Virus Preparation

### Safety Disclaimer

MICROSCOPY INNOVATIONS DOES NOT GUARANTEE THAT THIS PROTOCOL WILL RESULT IN THE COMPLETE DECONTAMINATION OF ALL VIRUSES. Follow all appropriate and necessary safety protocols for handling pathogenic viruses, including the United States Occupational Safety and Health Administration (OSHA) Bloodborne Pathogens Standard (29 CFR § 1910.1030). Users should also follow all appropriate safety guidelines for handling toxic reagents such as glutaraldehyde, uranyl acetate and osmium tetroxide.



**MICROSCOPY**  
INNOVATIONS<sub>LLC</sub>

213 Air Park Rd Suite 101  
Marshfield WI 54449-8626  
888.302.3925 or 715.384.3292  
[info@microscopyinnovations.com](mailto:info@microscopyinnovations.com)